



**Fowler
Simmons**
RADIOLOGY



InteleViewer™ InteleConnect™

Installation and Setup




InteleViewer™
by Intelerad

IntelViewer™ Software

The following steps explain how to install and setup IntelViewer™ software. If you already have IntelViewer™ installed, please proceed directly to Step 2.

Step 1: Install IntelViewer™

1. In your web browser, go to the following address:
<https://pacs.fsradiology.com.au/IntelViewer/>
(Please note this address is CASE sensitive)
2. Download the appropriate **IntelViewer™ Installer** (Windows 32-bit, Windows 64-bit or Mac OS X)
For Windows: Choose **Run** from the pop-up window to install.
For Mac: Double-click the downloaded file to install.
3. Accept the Terms and Conditions and click **Next**.
4. Click **Next** on any further screens that pop-up.
5. The **Install Complete** screen will appear when it has finished.




IntelViewer™ Installation Screen

Step 2: Setup IntelViewer™


1. Click the  IntelViewer™ icon to open the **Login Screen**.
2. In the **Login Screen**, click the  **Edit Bookmarks** button.
3. In the **Server Bookmarks** window, click the  **Add** button.
4. In the **New Bookmark** window, enter the following:
Server URL: <https://pacs.fsradiology.com.au>
(Do not hit ENTER / RETURN)
Ensure the **Protect Bookmark** tick-box is checked
Check the **Use Custom Server Name** tick-box
Name: Fowler Simmons
5. Click **Ok** to close the **New Bookmark** window.
6. Click **Close** to close the **Server Bookmarks** window.



IntelViewer™ Login Screen



IntelViewer™ Server Bookmarks Window



IntelViewer™ New Bookmark Window

IntelViewer™ Software

Step 3: Login to IntelViewer™

1. In the **Login Screen** enter the username and password provided to you by Fowler Simmons Radiology.
2. If you are presented with one of the following login errors, click the **IntelBrowser** link and proceed to Step 4.



IntelViewer™ Login Error 1

IntelViewer™ Login Error 2



IntelViewer™ Login Screen

Step 4: Change your Password

1. In the **IntelBrowser** screen, login with your username and current password.
3. Once you are logged-in, click **My Profile** in the menu and change your password.



IntelBrowser™ Window

IntelBrowser™ Menu

You are now ready to use the Fowler Simmons Radiology IntelViewer™ software.

If you would like to access images and reports without using the software, please see the following page for further information.

IntelConnect™ Online & App

Online Access

Access your images and reports through a web browser on any computer without installing the IntelViewer™ software.


1. In your web browser, go to the following address:
<https://pacs.fsradiology.com.au>
2. Login with your username and password.

iPhone or iPad App*

Access your images and reports on your mobile devices through the IntelConnect™ app on your iPhone or iPad.

1. Download the **IntelConnect™ app** from the Apple App Store on your iPhone or iPad. 
2. Open the app and enter the following details into the **New Account** screen:
Name: Fowler Simmons
Server: <https://pacs.fsradiology.com.au>
Username: *your username*
Password: *your password*
3. Click **Done** to add the account.
4. In the **Login Screen** enter your username and password.

*The IntelConnect™ app is currently only available for Apple devices.



IntelConnect™ Login Screen



Apple App Store

New Account Screen

Login Screen

Further Information and Enquiries

If you have not been provided with a username or password, or if you require help in setting up IntelViewer™, please contact:

Athanasia Savvoudis

Medical Liaison Officer
Phone 0429 855 722
asavvoudis@fsradiology.com.au

If you require assistance on how to use IntelViewer™, or have any clinical or Medical Imaging queries, please contact:

John Nisyrios

Chief Radiographer
Phone 08 8229 2100
jnisyrios@fsradiology.com.au